

## XSJ-2 Serials Microscope

Item	XSJ-2	XSX-1	XSX-2	XSY-1	Order NO.
Wide-field eyepiece WF10×-20mm	●●	●●	●●	●●	EW10
Wide-field eyepiece WF16×-14	●●	●●	●●	●●	EW16
WF10×-18 plan reticule eyepiece	●	●	●	●	EPr10b
8× centering eyepiece		●	●		EC8
4× plan achromatic objective, PL 4/0.10	●	●	●	●	B-OPr4
10×plan achromatic objective, PL(s) 10/0.25	●	●	●	●	B-OPr10
25×plan achromatic objective, PL(s) 25/0.40	●	●	●	●	B-OPr25
40×plan achromatic objective, PL(s) 40/0.65	●	●	●	●	B-OPr40
100×S-plan achromatic objective SPL(s) 100/1.25(Oil)	●	●	●	●	B-OSPor100
10×plan achromatic positive Ph-objective, PL(s) ph+ 10/0.25		●	●		B-OPpr10
40×plan achromatic positive Ph-objective, PL(s) ph+ 40/0.65		●	●		B-OPpr40
100×plan achromatic positive Ph-objective, SPL(s) ph+ 100/1.25(Oil)		●	●		B-OSPopr100
10×plan achromatic negative Ph-objective, PL(s) ph- 10/0.25			●		B-OPnr10
40×plan achromatic negative Ph-objective, PL(s) ph- 40/0.65			●		B-OPnr40
100×plan achromatic negative Ph-objective, SPL(s) ph- 100/1.25 (Oil)			●		B-OSPonr100
Main body, 6V30W	●	●	●	●	
Revolving trinocular head, 30° inclined	●	●	●	●	18X.30
Condenser N.A. 1.25	●	●	●	●	11X.10
Phase contrast Condenser		●	●		18X.2
135 photography device	●	●	●	●	11X.11
4×negative video eyepiece	●	●	●	●	EPHn4
Digital micro-photography device, for OLYMPUS C3000,C4000,C5000 serials	○	○	○	○	DV02
Green filterΦ38	●	●	●	●	11X.1.5-37
Blue filterΦ38	●	●	●	●	11X.1.5-38
Yellow filterΦ38	●	●	●	●	11X.1.5-39
Grey filterΦ38	●	●	●	●	11X.1.5-40
2.5×video eyepiece	○	○	○	○	EPHn2.5
6.3×video eyepiece	○	○	○	○	EPHn6.3
2.5×plan achromatic objective, PL 2.5/0.07	○	○	○	○	B-OPr2.5



Joy with your eyes

Chongqing MIC Technology Co., Ltd

Website: [www.micscope.com](http://www.micscope.com)

Email: [info@micscope.com](mailto:info@micscope.com)

TEL: +86-13436078184; FAX: +86-23-63913139

## XSJ-2 Serials Microscope

Dark filed condenser	○	○	○	○	18X.1
0.01mm micrometer 0.01mm	○	○	○	○	
Drawing device	○	○	○	○	18X.26
200W mercury lamp house				●	31X.2
200W mercury lamp power supply				●	31X.3
Epi-Flouresent illumination device, UV.V.B.G				●	18X.18
25×Flouresent objective FL 25/0.65				●	B-OFr25
40×Flouresent objective (oil) FL 40/1.00 GLYC				●	B-OFor40
100×Flouresent objective (oil) FL 100/1.25 GLYC				●	B-OFor100
C mount	○	○	○	○	CU43-2
Immersion oil				○	

Note: “●” in the table is standard attachment. “○” is optional accessories.

Design change: To keep pace with technological advances, we have reserved the right to make design modification and changes without notice.



Chongqing MIC Technology Co., Ltd  
 Website: [www.micscope.com](http://www.micscope.com)  
 Email: [info@micscope.com](mailto:info@micscope.com)

TEL: +86-13436078184; FAX: +86-23-63913139