

BA2000 Serials Microscope

Items	Spec and parameter	BA2202	BA2203	BA2302	BA2303	BA2500	Order No.
Eyepieces	10X-18	●●		●●			EP10a
	10X-20		●●		●●	●●	EW10
Plan Achromatic Objectives	SPlan10/0.25	●		●			B-OSPr10
	SPlan25/0.40	●		●			B-OSPr25
	SPlan40/0.65	●		●			B-OSPr40
	SPlan100/1.25(Oil)	●		●			B-OSPor100
	Plan10/0.25		●		●	●	B-OPr10
	Plan25/0.40		●		●		B-OPr25
	Plan40/0.65		●		●		B-OPr40
	Plan100/1.25(Oil)		●		●		B-OPor100
Epi-fluorescence Attachment	Fluorescence Objectives					●	B-OFr25
	25/0.65					●	B-OFr25
	40/1(Glycerin)					●	B-OFr40
	100/1.25(Glycerin)					●	B-OFr100
	Epi-fluorescence Main Body (B、G、V、UV)					●	FD06
	100W Mercury Lamp House					●	LH08
	Power Supply Case of 100W Mercury Lamp					●	LS14
	Glycerin					●	
Main Body	Quintuple Nosepiece	●	●	●	●	●	A08
Seidentopf Trinocular Head		●	●				TBR01
Seidentopf Trinocular Head				●	●	●	TTR01
Mechanical Stage	Moving Range: 76x50; Size: 160x135	●	●	●	●	●	S26
Condenser Focusing Team	Rack and Pinion Focusing Mechanism	●	●	●	●	●	CA03a
Condenser	N.A.1.25	●		●			C05
	N.A.1.25		●		●	●	C23
Field Diaphragm	Ø 2-ø35	●	●	●	●	●	DI06
Immersion Oil		●	●	●	●		
Filter	Ø 32(Blue)	●	●	●	●		CU02-1

Optional accessory

Item	Spec and Parameter	Order No.	
Eyepieces	10X-18 (Pointer)	Erg10	
	16X-10	EP16a	
	16X-14	EW16	
Plan Achromatic Objectives	S-Plan 4/0.10	B-OSPr4	
	Plan 4/0.10	B-Opr4	
Phase Contrast Attachment	Plan Positive Ph. Objective	S-Plan ph10/0.25	B-OSPpr10
		S-Plan ph40/0.65	B-OSPpr40a
		S-Plan ph100/1.25 (Oil)	B-OSPopr100
		Plan ph10/0.25	B-OPpr10
		Plan ph40/0.65	B-OPpr40
		Plan ph100/1.25/(Oil)	B-OPopr100
	Plan Negative Ph. Objective	S-Plan ph10/0.25	B-OPSnr10
		S-Plan ph40/0.65	B-OSPnr40a
		S-Plan ph100/1.25(Oil)	B-OSPonr100
		Plan ph10/0.25	B-OPnr10
		Plan ph40/0.65	B-OPnr40
		Plan ph100/1.25(Oil)	B-OPonr100
Ph. Condenser	Only for S-plan Ph. Objective	C06a	
	Only for Plan Ph. Objective	C06	
	Centering Eyepiece 8X	EC8	
Dark Field Attachment	Dark Field Condenser (N.A. 0.83-0.91)	C07	
Polarizing Attachment	Polarizer Slide	PO01	
	Analyzer Slide	PA01	
Photo Attachment	10X Reticule Eyepiece	EPr10b	
	2.5X Photo Eyepiece	EPHn2.5	
	4X Photo Eyepiece	EPHn4	
	6.3X Photo Eyepiece	EPHn6.3	
	135 Photo Device (With 4x Photo Eyepiece)	MP03	
TV and Video Attachment	C-mount Adaptor Color or Black and White Monitor CCD Camera	TV02	
Epi-illumination Attachment	6V 20W Halogen Tungsten Lamp (With Reflector)	II09	
Side Viewing Attachment	Side Viewing Attachment	TB06	
Epi-fluorescence Attachment	200W Mercury Lamp House	31X.2.2	
	Power Supply Case of 200W Mercury Lamp	31X.3	
Projection Attachment		PJ04	
Filter	Ø 32(Green、Amber)	CU02-2/3	
Image Analysis and Disposal System			



Chongqing MIC Technology Co., Ltd
Website: www.micscope.com
Email: info@micscope.com

TEL: +86-13436078184; FAX: +86-23-63913139

BA2000i Serials Microscope

Items	Spec and Parameter	BA2203i	BA2303i	BA2305i	Order No.
Eyepiece	WF 10X- Ø 20	● ●	● ●	● ●	EW10
Objectives	ISO Plan Achromatic Objectives 10X/0.25	●	●		B-OPi10a
	ISO Plan Achromatic Objectives 20X/0.40	●	●		B-OPi20a
	ISO Plan Achromatic Objectives 40X/0.66	●	●		B-OPi40a
	ISO Plan Achromatic Objectives 100X/1.25(Oil)	●	●		B-OPoi100a
	ISO Plan Fluorescent Objectives 10X/0.3			●	
	ISO Plan Fluorescent Objectives 25X/0.65			●	
	ISO Plan Fluorescent Objectives 40X/0.85			●	
	ISO Plan Fluorescent Objectives 100X/1.25 (oil)			●	
Main Body	Quintuple Nosepiece	●	●		A08a
Seidentopf Binocular Head		●			TBR11
Seidentopf Trinocular Head			●	●	TTR04
Mechanical Stage	Moving Range: 76×50 Size: 160×135	●	●		S26
Condenser Focusing Team	Rack and Pinion Focusing Mechanism	●	●		CA03a
Condenser	N.A. 1.25				C23
Field Diaphragm	Ø 2-Ø35	●	●		DI06
Immersion Oil		●	●		
Filter	Ø 32 (Blue)	●	●		CU2-1
Epi-fluorescent illumination system (V, UV, B, G), 100W AC mercury lamp house, 100W AC mercury power house (order number: FD06+LH08+PS-2)				●	

Optional accessory

Items	Spec and Parameter	Order No.	
Eyepieces	10X-18 (Pointer)	Erg10	
	16X-10	EP16a	
	16X-14	EW16	
IOS Plan Achromatic Objective	Plan 4X/0.10	B-OPi4a	
Phase Contrast Attachment	IOS Plan Positive Ph. Objectives	Plan ph10x/0.25	B-OPip10a
		Plan ph20x/0.40	B-OPip20a
		Plan ph40x/0.66	B-OPip40a
		Plan ph100x/1.25(Oil)	B-OPoip100a
	IOS Plan Negative Ph. Objectives	Plan ph10x/0.25	B-OPin10
		Plan ph20x/0.40	B-OPin20
		Plan ph40x/0.66	B-OPin40
		Plan ph100x/1.25(Oil)	B-Opoin100
	Phase Contrast Condenser		C22
	Centering Eyepiece		ECi
Dark Field Attachment	Dark Field Condenser (N.A. 0.83-0.91)	C07	
Polarizing Attachment	Polarizer Slide	PO01	
	Analyzer Slide	PA01	
Photo Attachment	10X Reticle Eyepiece	EPr10b	
	2.5X Photo Eyepiece	EPHn2.5	
	4X Photo Eyepiece	EPHn4	
	6.3X Photo Eyepiece	EPHn6.3	
	135 Photo Device (With 4x Photo Eyepiece)	MP03	
TV and Video Attachment	C-mount Adaptor、 Color or Black and White Monitor、 CCD Camera	TV02	
Epi-illumination Attachment	6V20W Halogen Tungsten Lamp (With Reflector)	H09	
Side Viewing Attachment	Two View Field Space 600mm	TB06	
Projection Attachment		PJ04	
Filter	Ø 32 (Green、 Amber)	CU02-2、 3	
Image Analysis and Disposal System			

Note: “●” in the table is standard attachment. “○” is optional accessories.

Design change: To keep pace with technological advances, we have reserved the right to make design modification and changes without notice.



Joy with your eyes

Chongqing MIC Technology Co., Ltd

Website: www.micscope.com

Email: info@micscope.com

TEL: +86-13436078184; FAX: +86-23-63913139