

BA3000i Microscope

Order NO.	Name	Item	BA3000i
EPh10b	10×High Point Wide Field Eyepiece/ WF10×-22		● ●
B-OPi4a	4×Infinity Plan Achromatic Objective/ PL 4×/-		●
B-OPi10a	10×Infinity Plan Achromatic Objective/ PL 10×/0.25		●
B-OPi20a	20×Infinity Plan Achromatic Objective/ PL(s) 20×/0.40		●
B-OPi40a	40×Infinity Plan Achromatic Objective/ PL(s) 40×/0.66		●
B-OPoi100a	100×Infinity Plan Achromatic Objective (Oil)		●
A15	Main Body/ Five Holes Transducer (two filling-in bits)		●
TTR07	Rotating Trinocular/ Binocular Adjustment (eyepiece tube sizeØ30)		●
S35	Stages/Fast Change/Size 180×142/Moving Range 80×55		●
CA11	Trays/Fast Change		●
C28	All light Achromatic Condenser/ N.A.0.92 /0.13		●
KI 01	Kola Illumination System/field irisØ2-Ø35		●
LH 09	Lamp House/12V 50W, Fuse 3A (Ø5x20)		●
11X.1.5-37	Green Filter/Ø38		●
11X.1.5-38	Blue Filter/Ø38		●
11X.1.5-39	Yellow Filter/Ø38		●
11X.1.5-40	Grey Filter/Ø38		●

Optional accessory

Order number	Name
EPhg10	Wide Field Eyepiece with Pointer (High Point) /WF 10×-22
EPh16b	High Point Wide Field Eyepiece /WF16×-16
B-OAi4b	Achromatic Objective 4×/-
B-OAi10b	Achromatic Objective 10×/0.25
B-OAi40b	Achromatic objective 40×/0.66 (s)
B-OAio100b	Achromatic Objective 100×/1.25 Oil (s)
B-OPip10a	Infinity Plan Positive Phase Objective /PL ph+ 10×/0.25
B-OPip20a	Infinity Plan Positive Phase Objective /PL(s) ph+ 20×/0.40
B-OPip40a	Infinity Plan Positive Phase Objective /PL(s) ph+ 40×/0.66
B-OPoip100a	Infinity Plan Positive Phase Objective /PL(s) ph+ 100×/1.25 Oil
B-OPin10	Infinity Plan Negative Phase Objective /PL ph- 10×/0.25
B-OPin20	Infinity Plan Negative Phase Objective /PL(s) ph- 20×/0.40
B-OPin40	Infinity Plan Negative Phase Objective /PL(s) ph- 40×/0.66
B-OPoin100	Infinity Plan Negative Phase Objective/PL(s) ph-100×/1.25 Oil



Joy with your eyes

Chongqing MIC Technology Co., Ltd

Website: www.micscope.com

Email: info@micscope.com

TEL: +86-13436078184; FAX: +86-23-63913139

ECia	Centering Telescope/20x/Ø30
C29	Phase Condenser /N.A.1.25
C30	Dark Field Condenser/ N.A. 0.83-0.91
PO05	Polarizer
PA04	Analyzer
EPhr10	High Point Wide Field Reticule Eyepiece 10×/Grid 0.1mm
EPHn2.5a	Photo Eyepiece 2.5×/ Image Size Ø42
EPHn4a	Photo Eyepiece 4×/ Image Size Ø42
MP06	135 Photo Device/ Match with 20 and 21
TV09b	Photo Device (With 0.35x System)/ 1/3"CCD
	0.01mm Micrometer/ Grid 0.01mm
TV09c+DV03	Digital Photo Device / OLYMPUS Digital Camera C5000 Serials
II 09	Reflective-transmitted Device/6V 20W
FD 06	Fluorescent Vertical Illumination /UV、 V 、 B、 G Four Kinds
B-OFi10	Infinity Plan Fluorescent Objective 10×/ FPL 10X/0.30
B-OFi25	Infinity Plan Fluorescent Objective 25×/ FPL 25X/0.65
B-OFi40	Infinity Plan Fluorescent Objective 40×/ FPL 40X/0.85
B-OFoi100	Infinity Plan Fluorescent Objective (glycerol) 100×/ FPL 100X/1.25 GLYC
LH08	DC Hg Lamp House 100W
PS-2	DC Hg Lamp House 100W
	Glycerol

Note: “●” in the table is standard attachment. “○” is optional accessories.

Design change: To keep pace with technological advances, we have reserved the right to make design modification and changes without notice.



Chongqing MIC Technology Co., Ltd
 Website: www.micscope.com
 Email: info@micscope.com

TEL: +86-13436078184; FAX: +86-23-63913139