

BK2000/3000 Serials Microscope

BK2000

Item	Specification	2100	2200	2201	2300	2301	Order No.
Eyepieces	Plan eyepiece (10X-18)	●	●●	●●	●●	●●	EP10a
	Photo negative eyepiece S4				●	●	EPHn4
Achromatic Objectives	4/0.10	●			○	○	B-OA4c
	10/0.25	●			○	○	B-OA10d
	(S) 40/0.65	●			○	○	B-OA40d
	(S) 100/1.25 oil	●			○	○	B-OAo100c
Plan objectives	S-Plan 4/0.10	○	●	●	●	●	B-OSPr4
	S-Plan 10/0.25	○	●	●	●	●	B-OSPr10
	S-Plan 40/0.65	○	●	●	●	●	B-OSPr40
	S-Plan (S) 100/1.25(oil)	○	●	●	●	●	B-OSPor100
Eyepiece Heads	Monocular head (45°inclined, 360°circumrotation)	●					TM18
	Teaching head	○	○	○			TM16
	Binocular head (45°inclined, slide)		●				TB10
	Trinocular head (45°inclined, slide)				●		TT05
	Binocular head (30°inclined, revolution)			●			TBR06
	Trinocular head (30°inclined, revolution)					●	TTR02
Main Body	Five-hole nosepiece (built-in)		●	●	●	●	MSBk4
	Four-hole nosepiece (built-in)	●	○	○	○	○	MSBk3
Condenser	N.A.1.25Oil	●	●	●	●	●	C20
Base	AC220V or 110V	●	●	●	●	●	B26
Photo Attachment	35mm DF camera				●	●	
	Photo tube				●	●	MP03
Immersion oil		●	●	●	●	●	
Filter	Ø32 (blue)	●	●	●	●	●	CU02-1
CCD attachment					○	○	TV02

Optional accessory

Item	Specification	Order No
Eyepieces	Plan eyepiece (16X-10)	EP16a
	Pointer eyepiece (10X-17.2)	Erg10
	Plan reticle eyepiece (10X-18)	EPr10b
	2.5X photo negative eyepiece	EPHn2.5
Achromatic Phase-contrast Attachment	8X Central eyepiece	EC8
	Achromatic negative phase-contrast objective 10/0.25	B-OAn10
	Achromatic negative phase-contrast objective 40/0.65	B-OAn40
	10x annular plate	RP01
	40x annular plate	RP02
Plan Phase-contrast Attachment	Plan achromatic negative phase-contrast objective SPlan10/0.10	B-OSPnr10a
	Plan achromatic negative phase-contrast objective S-Plan (S) 40/0.10	B-OSPnr40a
	Plan achromatic negative phase-contrast objective S-Plan (S) 100/1.25	B-OSPnr100
	Plan achromatic positive phase-contrast objective SPlan10/0.10	B-OSPpr10
	Plan achromatic positive phase-contrast objective S-Plan (S) 40/0.65	B-OSPpr40
	Plan achromatic positive phase-contrast objective S-Plan (S) 100/1.25	B-OSPopr100
	Phase-contrast condenser (used in S-Plan phase-contrast Objective)	C06a
Field diaphragm		DI06
Filter	Ø 32 (Green)	CU02-2
	Ø 32 (Amber)	CU02-3
Reflector		MT04
Dark field Attachment	Dark field condenser N.A.0.83-0.91	C07
	Simple dark field attachment	DS03
Simple polarizer Attachment	Polarizer slide	PO01
	Analyzer	PA01

BK3000

Item	Specification	BK3200	BK3300	Order No.
Eyepieces	Plan eyepiece (10X-18)	● ●	● ●	EP10a
	Photo negative eyepiece S4		●	EPHn4
IOS achromatic Objectives	4X/0.10	●	●	B-OAi4a
	10X/0.25	●	●	B-OAi10a
	(S) 40X/0.66	●	●	B-OAi40a
	(S) 100X/1.25(oil)	●	●	B-OAoil100a
Eyepiece Heads	Binocular head(30°inclined,revolution)	●		TBR11
	Trinocular head(30°inclined,revolution)		●	TTR04
Main body	Five-hole nosepiece	●	●	MSBk5
Condenser	N.A.1.25Oil	●	●	C20
Base	AC220V or 110V	●	●	B26
Photo Attachment	Photo tube		●	MP03
	35mm DF camera		●	
Immersion oil		●	●	
Filter	Ø 32 (blue)	●	●	CU02-1

Optional accessory

Item	Specification	Order No.
Eyepieces	Plan eyepiece (16X-10)	EP10a
	Pointer eyepiece (10X-17.2)	Erg10
	Plan reticle eyepiece (10X-18)	EPr10b
	2.5X photo negative eyepiece	EPHn2.5
Filter	Ø 32 (Green)	CU02-2
	Ø 32 (Amber)	Cu02-3
Reflector		MT04
Dark field Attachment	Dark field condenser	C07
	Simple dark field attachment	DS03
Simple polarizer Attachment	Polarizer slide	PO01
	Analyzer slide	PA01

Note: “●” in the table is standard attachment. “○” is optional accessories.

Design change: To keep pace with technological advances, we have reserved the right to make design modification and changes without notice.



Joy with your eyes

Chongqing MIC Technology Co., Ltd

Website: www.micscope.com

Email: info@micscope.com

TEL: +86-13436078184; FAX: +86-23-63913139