

BL Serials Microscope

Model		BL 101	BL 102	BL 103	BL 104	BL 201	BL 201	Order No
Items Spec. and parameter								
Eyepieces	Ramsden 10X	●	●	●	●	●	●	ER10
	Pointer 10X	○	○	○	○	○	○	EP10f
Objectives	Achromatic 4X/0.10	●	●	●	●	●	●	B-OA4b
	Achromatic 10 X /0.25	●	●	●	●	●	●	B-OA10b
	Achromatic 40 X /0.65	●	●	●	●	●	●	B-OA40b
	Achromatic 100 X /1.25 (oil)		○		○		○	B-OA100 b
Head	Monocular head (inclined 45°)	●	●					130X.1
	Teaching head (inclined 45°)			●	●			130X.5
	Binocular head (inclined 45°)					●	●	130X.6
Body	Standard main body	●	●			●	●	130X.2
	Teaching main body			●	●			130X.2A
Nosepiece	Build in quadruple nosepiece	●	●	●	●	●	●	
	Build in quindruple nosepiece							
Stage	Fixed Stage	●		●		●		130X.3
	Mechanical Stage		●		●		●	
Condenser	Single condenser, 5 aperture disk diaphragm	●		●		●		
	ABBE condenser		●		●		●	
Base	220V50W fluorescent lamp	●	●	●	●	●	●	130X.4
LED light resource		○	○	○	○	○	○	
Stage clip		●●		●		●		CL 05

Note: “●” in the table is standard attachment. “○” is optional accessories.

Design change: To keep pace with technological advances, we have reserved the right to make design modification and changes without notice.



Chongqing MIC Technology Co., Ltd
 Website: www.micscope.com
 Email: info@micscope.com

TEL: +86-13436078184; FAX: +86-23-63913139