

BP Serials Microscope

Model		BP 101	BP 102	BP 103	BP 104	BP 105	BP 201	Order No
Items Spec. and parameter								
Eyepieces	Ramsden WF10X /ø17.2	●	●	●	●	●	●	ER10
	Ramsden (pointer) WF10X /ø17.2	○	○	○	○	○	○	ERg10
	Huygenian 10X	○	○	○	○	○	○	EH 10
Objectives	Achromatic 4X/0.10	●	●	●	●	●	●	8604
	Achromatic 10X/0.25	●	●	●	●	●	●	8601
	Achromatic (S) 40X/0.65	●	●	●	●	●	●	8209A
	Achromatic (S) 100X/1.25(OIL)				●		●	8210A
Head	Monocular head (inclined 45°)	●	●	●	●			Bp05
	Monocular head (inclined 45°)					●		Bp05a
	Binocular head (inclined 45°)						●	X1A
	Teaching head	○	○	○	○		○	Bp09
Body	Only coarse focusing attachment	●						Bp003
	Only coarse focusing attachment					●		Bp003a
	Separate coarse and fine focusing attachment		●	●	●		●	Bp006
Stage	Match Abbe condenser with leaver				●		●	S25b
	Single condenser, 5 aperture disk diaphragm			●		●		S02a
	Single condenser, match iris diaphragm	●						Bp07a
	Single condenser, 5 aperture disk diaphragm		●					S25
Base	Mirror (ø45mm)					●		B16c
	Incandescent lamp (200V, 20W)		●	●		●		B16d
	Fluorescent lamp (110V, 5W)	●			●			B16e
Exterior fixing Triple nosepiece		●	●	●				N/R3
Exterior fixing Quadruple nosepiece					●	●	●	N/R4
Iris Diaphragm		●		●				8411
Condenser					●		●	C21
Stage clip		●●	●●	●●	●●	●●		CL05
Attachable mechanical stage							●	Ed01
Hexagon spanner (M2)		●	●	●	●	●	●	

Note: “●” in the table is standard attachment. “○” is optional accessories.

Design change: To keep pace with technological advances, we have reserved the right to make design modification and changes without notice.



Joy with your eyes

Chongqing MIC Technology Co., Ltd

Website: www.micscope.com

Email: info@micscope.com

TEL: +86-13436078184; FAX: +86-23-63913139