

T-1 Serials Microscope

Name	Specifications	T-1	T-1LED
Viewing head	Monocular head, 45°inclined, 360°rotating	●	●
Eyepiece	Wide-field eyepiece WF10×-18mm	●	●
	Wide-field eyepiece WF10×-18mm with pointer	○	○
	Wide-field eyepiece WF10× with scale	○	○
	Wide-field eyepiece WF16×	○	○
	Wide-field eyepiece WF20×	○	○
Nosepiece	Triple nosepiece with positive click stopper	●	●
	Quadruple nosepiece with built-in click stopper	○	○
Objective	DIN achromatic objective 4× / NA0.10	●	●
	DIN achromatic objective 10× / NA0.25	●	●
	DIN achromatic objective 20×R / NA0.40 (spring loaded)	○	○
	DIN achromatic objective 40×R / NA0.65 (spring loaded)	●	●
	DIN achromatic objective 60×R / NA0.85 (spring loaded)	○	○
	DIN achromatic objective 100×R(oil) / NA1.25 (spring loaded)	○	○
Stage	Mechanical stage 110mm×120mm with paired clips	●	○
	Mechanical stage 110mm×120mm with movable specimen holder	○	●
Condenser	N.A.0.65 condenser, disc diaphragm	●	●
Illuminati-on	An incandescent bulb 20W/115V or 20W/230V.	●	○
	LED lamps, Ni-Cd rechargeable battery pack, an AC/DC adaptor. The brightness is adjustable.	○	●
	Halogen bulb 10W/12V and a power adaptor for 115V or 230V. The brightness is adjustable.	○	○
Focusing system	Coaxial focusing system, with safe rack stopper	●	●
Optional Accessories	Movable specimen holder (black or slive)	○	○
	DCM35 (350K pixels, USB2.0) digital camera for microscope	○	○
	DCM130 (1.3M pixels, USB2.0) digital camera for microscope	○	○
	Wooden case with lock & keys	○	○
	Aluminium-Alloy Enstrument Case	○	○
Packing	Each set in a foam in an inner carton, 4 inner cartons in a master carton.	●	●

Note: “●” in the table is standard attachment. “○” is optional accessories.

Design change: To keep pace with technological advances, we have reserved the right to make design modification and changes without notice.

**CE ISO9001
ISO14001**

MIC[®]

Joy with your eyes

Chongqing MIC Technology Co., Ltd

Website: www.micscope.com

Email: info@micscope.com

TEL: +86-13436078184; FAX: +86-23-63913139