

## XDS-1C Inverted Biological Microscope

Order No.	Items	XDS-1C
EW10	Eyepiece WF10×/ø20	●●
EW16	Eyepiece WF16×/ø14	●●
EPr10b	Plan10×/ø18 reticle Eyepiece	●
B-OLPr10	LWD Plan Achromatic Objective 10/0.25,160/1.5,W.D.=9.36mm	●
B-OLPr25	LWD Plan Achromatic Objective 25/0.40,160/1.5,W.D.=6.50mm	●
B-OLPr40	LWD Plan Achromatic Objective 40/0.65,160/1.5,W.D.=4.50mm	●
B-OLPpr10	LWD Plan Positive Phase Contrast Objective 10/0.25,160/1.5,W.D.=9.36mm	○
B-OLPpr25	LWD Plan Positive Phase Contrast Objective 10/0.40,160/1.5,W.D.=6.50mm	●
B-OLPpr40	LWD Plan Positive Phase Contrast Objective 40/0.65,160/1.5,W.D.=4.50mm	○
B-OLPnr10	LWD Plan Negative Phase Contrast Objective 10/0.25,160/1.5,W.D.=9.36mm	○
B-OLPnr25	LWD Plan Negative Phase Contrast Objective 10/0.40,160/1.5,W.D.=6.50mm	○
B-OLPnr40	LWD Plan Negative Phase Contrast Objective 40/0.65,160/1.5,W.D.=4.50mm	○
TBR03	Seidentopf Binocular Head (inclined 45°) Interpulary distance: 55~75mm	●
S09A	Mechanical Stage Size: 180×155mm; Moving range:75×50mm	●
33X-4	Set for Stage Insert Plate	●
33X-5	Stage Insert plate (V)	●
9J-8	Stage Insert plate (III)	○
9J-9	Stage Insert plate (IV)	●
RP09	Annular Spot 10×	○
RP10	Annular Spot 25×	●
RP11	Annular Spot 40×	○
9J.7	Photo Adaptor With Photo Eyepiece S5×	○
CL02	Stage Clips	●●
EPHp5	5× video eyepiece	○
EPHp6.3	6.3×video eyepiece	○
9J.7B	Photograpy Device	○
9J.15+Cu43-1+	Photograpy adaptor      Photograpy Device	○
33X.3c	Fixed Illumination system	●
CU02-4~7	Blue, Green Yellow, Gray filterΦ32	●
	Micrometer (division 0.01mm)	○

Note: “●” in the table is standard attachment. “○” is optional accessories.

Design change: To keep pace with technological advances, we have reserved the right to make design modification and changes without notice.



Joy with your eyes

Chongqing MIC Technology Co., Ltd

Website: [www.micscope.com](http://www.micscope.com)

Email: [info@micscope.com](mailto:info@micscope.com)

TEL: +86-13436078184; FAX: +86-23-63913139