

XSP-M Serials Microscope

Model		M1	M2	M16	M2A	Order No
Items Spec. and parameter						
Eyepieces	Ramsden 10X	●	●	●	●	ER10
	Ramsden 10X (Pointer)	○	○	○	○	ERg10
Objectives	Achromatic 4X/0.10	●	●	●	●	B-OA4a
	Achromatic 10 X /0.25	●	●	●	●	B-OA10a
	Achromatic 40 X /0.65 (spring)	●	●	●	●	B-OA40a
	Achromatic 100 X /1.25 (spring, oil)	○	○	●	●	B-OA100a
Head	Monocular head (inclined 45°)	●	●	●	●	TM10
	Teaching head	○	○	○	○	TM05a
Body	Standard main body	●	●	●	●	MS01
Nosepiece	Triple nosepiece	●	●			N/R4
	Quadruple nosepiece			●	●	
Stage		●	●	●	●	
Condenser				●	●	2XB-14
Moving scale		○	○	●	●	
Base		●	●	●	●	110X0.1
Clip		●●	●●		●●	CL05

Note: “●” in the table is standard attachment. “○” is optional accessories.

Design change: To keep pace with technological advances, we have reserved the right to make design modification and changes without notice.



Chongqing MIC Technology Co., Ltd
 Website: www.micscope.com
 Email: info@micscope.com

TEL: +86-13436078184; FAX: +86-23-63913139