

## XSZ-D2 Microscope

Order NO.	Items	XSZ-D2
EW10	Plan Eyepiece WF10 / 20mm	● ●
EW16	Plan Eyepiece WF16 / 14mm	● ●
EPr10c	Plan Reticle Eyepiece PF10× / 18mm	●
B—OLPr10	LWD Plan Achromatic Objective 10× / 0.25 160 / 1.5 WD.=7.9mm	●
B—OLPr25	LWD Plan Achromatic Objective 25× / 0.40 160 / 1.5 WD.=5mm	●
B—OLPr40	LWD Plan Achromatic Objective 40× / 0.65 160 / 1.5 WD.=3mm	●
B—OLPnr25	LWD Plan Positive Phase Objective 25× / 0.40 160 / 1.5 WD.=5mm	●
EC5	Centering eyepiece	●
12X.2.3-19,20	25X Phase Contrast Annular	●
	Sliding Binocular Head (Inclined 45°)	●
12X-3	Stage (On The Main Body) Size:170X166, moving range: 20X20	
8204	Stage (On The Main Body) Size:175X160, moving range: 80X50	●
	Videoing group	●
EPHn.4	4X video negative eyepiece	●
	Halogen Tungsten Lamp 6V / 20W	● ● ●
7X.5.1	Cultural Vessel Bottom Thickness 1.5mm	●
7X.5.2	Cultural Flask (Bottom Thickness 1.5mm)	●

Note: “●” in the table is standard attachment. “○” is optional accessories.

Design change: To keep pace with technological advances, we have reserved the right to make design modification and changes without notice.



Joy with your eyes

Chongqing MIC Technology Co., Ltd

Website: [www.micscope.com](http://www.micscope.com)

Email: [info@micscope.com](mailto:info@micscope.com)

TEL: +86-13436078184; FAX: +86-23-63913139