

## MA3000 Serials Microscope

Order No.	Item	MA3000
15Q	MA3000 Reflecting Main Body For Industrial Inspection Reflecting, Manual Nosepiece, Quick Change Set Device Of Bright Field And Dark Field, Polarizing Attachment	●
S37	200mm Moving Stage Size: 375mmX335mm, Moving Range: 200mmX200mm; With Quick Cluch	●
15Q-4	Glass Stage 200mmX200mm	●
LH01	Lamp Housing / Reflected Light Source (12V50W)	●
TTR07a	Seidentopf Trinocular Inner Dia. of The Eyepiece Tube is $\phi 30$ , Inclined 30°	●
	Quintuple Manual Nosepiece Two Big Holes for Dark Field Objective, Three Small Holes For Bright Field Objective	●
M-OLPi4	4X Long Distance Plan Objective Infinite aberration correction WD=27.236mm NA=0.1	●
M-OLPi10	10X Long Distance Plan Objective Infinite aberration correction WD=18.48mm NA=0.25	●
M-OLPi20	20X Long Distance Plan Objective Infinite aberration correction WD=8.35 mm NA=0.4	●
M-OLPi40	40X Long Distance Plan Objective Infinite aberration correction WD=3.895mm NA=0.65	●
M-OLPi50	50X Long Distance Plan Objective Infinite aberration correction WD=1.95mm NA=0.7	●
M-OLPi80	80X Long Distance Plan Objective Infinite aberration correction WD=0.848mm NA=0.8	●
M-OLPid10	10X Long Distance Plan Dark Field Objective Infinite aberration correction WD=18.48mm NA=0.25	●
M-OLPid20	20X Long Distance Plan Dark Field Objective Infinite aberration correction WD=8.35mm NA=0.4	●
M-OLPid40	40X Long Distance Plan Dark Field Objective Infinite aberration correction WD=3.895mm NA=0.65	●
	50X Long Distance Plan Dark Field Objective Infinite aberration correction WD=1.95mm NA=0.7	
	80X Long Distance Plan Dark Field Objective Infinite aberration correction WD=0.848mm NA=0.8	
Eph10b	10X Wide Field Objective FN= $\phi$ 22mm	●●
15Q.1.6-9	Yellow, Green Blue, Neutral Filter (One Piece For Each)	●
15Q.1-10	Ground Glass	●
EPHn2.5	2.5X Negative Photo Eyepiece for SLR(Single Lens Reflex) camera only	●
MP05	Digital Adapter C Mount Suit to DV02 or Video Camera	●
DV02	Can Be Used with OLYMPUS5000 Serials Camera	●
	Power Line	●



Joy with your eyes

Chongqing MIC Technology Co., Ltd

Website: [www.micscope.com](http://www.micscope.com)

Email: [info@micscope.com](mailto:info@micscope.com)

TEL: +86-13436078184; FAX: +86-23-63913139

## MA3000 Serials Microscope

5X20	Fuse	●
EW16b	16X Wide Field Objective FN=ø 16mm	○○
	Semiconductor crystal clip, ø200mm semiconductor crystal can be placed on it	○
MP06	35mm video tube	○
EPHn4	Photo eyepiece 4X	○
	Photo eyepiece 6.3X	○

Note: “●” in the table is standard attachment. “○” is optional accessories.

Design change: To keep pace with technological advances, we have reserved the right to make design modification and changes without notice.



Chongqing MIC Technology Co., Ltd  
Website: [www.micscope.com](http://www.micscope.com)  
Email: [info@micscope.com](mailto:info@micscope.com)

TEL: +86-13436078184; FAX: +86-23-63913139