

IBE-2000 Microscope

Order	Item		Configuration		
			IBE2000	IBE2003	Order No
1	Wide field eyepiece	10× F.N=20	●●	●●	EW10
2		16× F.N=14	●●	●●	EW16
3	Long working distance plan achromatic objectives	2.5×	○	○	B-OLPr2.5
4		4×	○	○	B-OLPr4
5		10× WD=7.9	●	●	B-OLPr10
6		25× WD=5	●	●	B-OLPr25
7		40× WD=3	●	●	B-OLPr40
8	Long working distance plan achromatic positive phase contrast objectives	10× WD=7.9	●	●	B-OLPpr10
9		25× WD=5	●	●	B-OLPpr25
10		40× WD=3	●	●	B-OLPpr40
11	Long working distance plan achromatic negative phase contrast objectives	10× WD=7.9	○	○	B-OLPnr10
12		25× WD=5	○	○	B-OLPnr25
13		40× WD=3	○	○	B-OLPnr40
14	Long working distance phase contrast condenser	Mode I ,W.D=21.4, NA0.55	○	○	C13
15		Mode II, W.D=32, NA0.4	○	○	C14
16		Mode III, W.D=55, NA0.3	●	●	C15
17	Centering eyepiece	20×centering eyepiece	●	●	Eci20
18	Reticule eyepiece	10×plan reticule eyepiece	●	●	Epr10b
19	Body	IBE2000 main body	●	●	
20	Eyepiece unit	Binocular (Trinocular) unit	●	●	18X.30
21	Photographic equipment	Photographic adapter	●	●	MP04
22		35mm SLR camera, shutter release	●	●	
23		Digital camera adapter	○	○	DV02
24	CCTV equipment	Video camera adapter	○	○	TV04
25		“C”mount adapter	○	○	CU43-2
26	Culture bottle and vessel	Culture bottle and vessel	●	●	7X.5.1/7X.5.2
27	Incidence fluorescence equipment	EPI-fluorescence illuminator(UV, V, B, G)		●	FD05
28		EPI-fluorescence illuminator (B,G)	●	●	FD05a
29		Centering objective	●	●	CU61
30		100W lamp house	●	●	LH08
31		DC 100W power supply unit	●	●	LS14
32	Micrometer	0.01mm section	●	●	

Note: “●” in the table is standard attachment. “○” is optional accessories.

Design change: To keep pace with technological advances, we have reserved the right to make design modification and changes without notice.



Chongqing MIC Technology Co., Ltd
 Website: www.micscope.com
 Email: info@micscope.com

TEL: +86-13436078184; FAX: +86-23-63913139