

L-I1 Serials Microscope

Model	L-I1
Eyepiece	Wide field WF10X(μ 20mm)
Objective	Long working distance plan achromatic PL L 10X/0.25
	Long working distance plan achromatic PL L 25X/0.40
	Long working distance plan achromatic PL L 40X/0.60
	Long working distance plan achromatic phase contrast objective PL L 10X/0.25 PHP
Eyepiece tube	Trinocular (Inclination of 30°, 100% image light for photography capable)
Focus system	Coaxial coarse/fine focus, with tensional adjustable and up stop minimum division of fine focusing: 2 μ m.
Nosepiece	Quadruple(Ball bearing inner locating)
Stage	Double layer mechanical(Size200mmX152mm,moving range:75mmX30mm)
Shelf of specimen	Match to install μ 68mm or 77mmX33mm culture utensil
	Match to install μ 68mm culture utensil
Filter	Blue filter
	Green filter
	Frosted filter
Impellent style phase contrast unit	Centering telescope
	10X Impellent style phase plate
	Condenser Rack & pinion adjustable, NA=0.4,working distance: 30mm
Light source	6V 20W Halogen lamp with brightness control.



Joy with your eyes

Chongqing MIC Technology Co., Ltd

Website: www.micscope.com

Email: info@micscope.com

TEL: +86-13436078184; FAX: +86-23-63913139

L-I1 Serials Microscope

Optional accessories

Name	Sort/Technique parameter	NO.	
Eyepiece	WF16X(μ 11mm)	1051016	
	Dividing 10X(μ 20mm), 0.1mm/Div	1120010	
Nosepiece	Quintuple(Ball bearing inner locating)	032001	
Condenser	Ultra long working distance, W.D.:70mm	050006	
Filter	Yellow filter	115003	
Impellent style phase contrast unit	Objectives	Long working distance plan achromatic phase contrast objective PL L 25X/0.40 PHP	2430125
		Long working distance plan achromatic phase contrast objective PL L 40X/0.60 PHP	2430140
	Impellent phase plate	25X	112002
		40X	112003
Sliding phase contrast unit	Objectives	Long working distance plan achromatic phase contrast objective PL L 10X/0.25 PHP2	2430310
		Long working distance plan achromatic phase contrast objective PL L 25X/0.40 PHP2	2430325
		Long working distance plan achromatic phase contrast objective PL L 40X/0.60 PHP2	2430340
	Sliding phase plate	10X	112004
		25X/40X	112005
	Condenser Rack & pinion adjustable, NA=0.4, working distance: 55mm		050007
CCD adapter	0.4X	812001	
	0.5X	812004	
	1X	812002	
	0.5X with dividing 0.1mm/Div	812003	
Camera	DV-1(With USB & video output)	800001	
	DV-2(With USB output)	800003	
	DV-3(With video output)	800005	
Digital camera adapter	CANON digital camera adapter(A610,A620,A630,A640)	821001	
Photograph unit	2.5X/4X Change over photograph attachment with 10X viewing eyepiece	840001	
	4X Focusing photograph attachment	840002	
	MD Adapter	840003	
	PK Adapter	840004	

Note: “●” in the table is standard attachment. “○” is optional accessories.

Design change: To keep pace with technological advances, we have reserved the right to make design modification and changes without notice.



Chongqing MIC Technology Co., Ltd
 Website: www.micscope.com
 Email: info@micscope.com
 TEL: +86-13436078184; FAX: +86-23-63913139