

## NY-1 Serials Microscope

Item	Specification		
Optical system	Infinite Optical System	●	
Main body	NY-1	●	
Head	Seidentopf binocular Head (Inclined 45°)	●	
Eyepiece	10X/Φ22mm Wide field eyepiece	●	
	15X/Φ16mm Wide field eyepiece	○	
	20X/Φ12mm Wide field eyepiece	○	
	10X/Φ20mm Reticle eyepiece(0.1mm)	●	
Nosepiece	Quintuple Nosepiece	●	
Objective	Infinite long workingdistance plan achromatic objective (bright and dark field )	5X/0.12 WD: 10mm	●
		10X/0.25 WD: 10mm	●
		20X/0.4 WD: 5mm	●
		50X/0.75 WD: 1.3mm	●
		100X/0.90 WD: 0.7mm	●
Focusing system	Coaxial focusing: 0.002mm	●	
	Moving range(from the focus of the stage) 1mm upside,7mm downside	●	
Stage	Size: 226mm×178mm	●	
	Slide holder two piece, Small stage:	●	
Mechanical Stage	Moving range: 40mm×40mm	●	
Illumination	12V50 Halogen Lamp	●	
DIC	Applicable with 5X、 10X、 20X、 50X、 100X objective	○	
Filter	Inner turntable (green, blue, yellow, gray, gross glass)	●	
Accessory	Polarizer, micrometer(0.01mm), Centring objective	●	
	Digital camera adaptor	○	
	CCD adaptor	○	
	135 video objective	○	

Note: “●” in the table is standard attachment. “○” is optional accessories.

Design change: To keep pace with technological advances, we have reserved the right to make design modification and changes without notice.



Joy with your eyes

Chongqing MIC Technology Co., Ltd

Website: [www.micscope.com](http://www.micscope.com)

Email: [info@micscope.com](mailto:info@micscope.com)

TEL: +86-13436078184; FAX: +86-23-63913139