

D-M1 Microscope

		D-M1
Viewing Head	Trinocular head inclined at 30°, interpupillary distance 48-75mm	●
Eyepiece	High-point , extra wide field eyepiece EW10×/ 20	●
Infinite Plan Achromatic Objective	4×/ 0.1	●
	10×/ 0.25	●
	20×/ 0.40	●
	40×/ 0.60	●
	50×/ 0.75	○
	80×/ 0.90	○
Nosepiece	Quintuple nosepiece	●
Stage	Plain stage	●
	Attachable mechanical stage, X-Y coaxial control, moving range 120×78mm	○
	Auxiliary stage	○
Focusing	Coaxial coarse & fine adjustment, Vertical objective movement, Coarse stroke 37.7mm per rotation, Fine stroke 0.2mm per rotation	●
Kohler Illumination	Halogen lamp 6V/ 30W	●
Filter	Blue, Yellow, Green and Ground Glass	●
Polarization Attachment		●
Photo Attachment		○
Video Adapter	Video adapter with C mount	○

Note: “●” in the table is standard attachment. “○” is optional accessories.

Design change: To keep pace with technological advances, we have reserved the right to make design modification and changes without notice.



Chongqing MIC Technology Co., Ltd
 Website: www.micscope.com
 Email: info@micscope.com

TEL: +86-13436078184; FAX: +86-23-63913139