

H-ZM2 Microscope

Objectives:	4*0.1W.D.22mm,10*0.25W.D.10mm,20*0.4W.D.5mm, 40*0.6W.D.3mm
Eyepiece	WF10X/22mm
Viewing head	compensation free binocular head
Nosepiece	Quintuple nosepiece
Stage	Double layers mechanical stage :190*160mm, moving stage:75*50mm
Condenser	NA1.25Abbe condenser with iris diaphragm
Focusing	Coaxial coarse&fine focusing adjustment with rack and pinion mechanism fine focusing scale value 0.002mm
Illumination	Epi-kohler illumination.with aperture iris diaphragm and field iris diaphragm transmission light,power
Power	220V/110V Halogen lamp,Adjustable brightness
Filter	green:blue:yellow.
Parts for choice	Polarizing,Polarizing outfit

Note: “●” in the table is standard attachment. “○” is optional accessories.

Design change: To keep pace with technological advances, we have reserved the right to make design modification and changes without notice.

CE **ISO9001**
ISO14001



Chongqing MIC Technology Co., Ltd
Website: www.micscope.com
Email: info@micscope.com

TEL: +86-13436078184; FAX: +86-23-63913139