

MX1000 Microscope

Order NO.	Item	MX1000
EPh10b	10× high eye point wild-field eyepiece WF10×/22 f=25 eye-pointdistance 25	● ●
EPh16b	16×high eye point wild-field eyepiece WF16×/16 f=15.6 eye-pointdistance 18	● ●
EPhr10	10×high eye point wild-field reticle eyepiece WF10×/22 f=25	●
M-OLPi10	10×infinite long distance plan objective NA=0.25 WD=18.48	●
M-OLPi20	20×infinite long distance plan objective NA=0.40 WD=8.35	●
M-OLPi40	40×infinite long distance plan objective NA=0.65 WD=3.895	●
M-OLPi50	50×infinite long distance plan objective NA=0.70 WD=1.95	●
M-OLPi80	80×infinite long distance plan objective NA=0.80 WD=0.848	●
M-OLPid10	10×infinite long distance plan dark field objective NA=0.25 WD=18.48	●
M-OLPid20	20×infinite long distance plan dark field objective NA=0.40 WD=8.35	○
M-OLPid40	40×infinite long distance plan dark field objective NA=0.65 WD=3.895	●
M-OLPid50	50×infinite long distance plan dark field objective NA=0.70 WD=1.95	○
M-OLPid80	80×infinite long distance plan dark field objective NA=0.80 WD=0.848	○
S28	Stage---Size: 238×160 Moving range: 55×52	●
S28.1	Slide: external circleφ110	●
TTR22	Seidentopf trinocular head 30° inclined	●
PO02	Polarizer	●
PA01a	Analyzer	●
CL02	Sample presser group	● ●
TV09b	Digital video device	○
TV04a	Digital video device	●
SC03	Sample clip	○
LH15	Lamp house	●
	Tungsten 24V/100W	● ●
	LCD liquid crystal display	○
	VGA digital camera	○
	0.01mm micrometer	○

MX1000 Microscope

Note:

“●” in the table is standard attachment. “○” is optional accessories.

Design change: To keep pace with technological advances, we have reserved the right to make design modification and changes without notice.



www.micscope.com



Chongqing MIC Technology Co., Ltd

Website: www.micscope.com

Email: info@micscope.com

TEL: +86-13436078184; FAX: +86-23-63913139