

NY-2 Serials Microscope

Item	Specification	NY-2
Main body		●
Objective	PC4X/0.10	●
	PC10X/0.20	●
	PC25X/0.40	●
	PC40X/0.65	●
	PC63X/0.85	●
Eyepiece	10X	●
	10X (0.1mm reticule Eyepiece)	●
magnification system	0.8X	●
	1X	●
	1.25X	●
Video eyepiece	4X	●
Dark field collar	PC10X/0.20	●
	PC25X/0.40	●
Glass field iris diaphragm	φ13mm	●
	φ21mm	●
	φ27mm	●
Metal circular orifice iris diaphragm	φ10mm	●
	φ15mm	●
	φ20mm	●
Objective micrometer	0.01mm	●
Automatic shutter control		●
cassette (plate-holder)		●
Imaging plate	Gross glass with centre line	●
	Gross glass with 100mm reticle	●
	Gross glass with X reticle	●
	Gross glass with 8 degree grain size	●

NY-2 Serials Microscope

Focusing glass	6X	●
Polarizer		●
Analyzer		●
Light	12V50W Halogen Lamp	●
Light screen		●
135 photography device		○
CTV video adaptor		○
Cable		●
Cable release		●
Stage tableting		●
Optional	Eyepiece 12.5X、 16X, digital camera adaptor, image analysing systme, DF135 camera, digital video accessory	

Note: “●” in the table is standard attachment. “○” is optional accessories.

Design change: To keep pace with technological advances, we have reserved the right to make design modification and changes without notice.



Chongqing MIC Technology Co., Ltd
 Website: www.micscope.com
 Email: info@micscope.com

TEL: +86-13436078184; FAX: +86-23-63913139