

## O-I1 Inverted Metallurgical Microscope

Items	Specification	O-I1
Eyepieces	WF10×/Φ18	●●
	WF10×/Φ20	○○
	WF12.5×/Φ14	●●
	WF16×/Φ13	○○
	Huygenian Ocular 5×/Φ20	○○
	WF10×/Φ18 (Reticule)	●
Plan Achromatic Objectives	Plan 2.5×/0.07	○
	Plan 4×/0.10	○
	Plan 10×/0.25	●
	Plan 20×/0.35(S)	●
	Plan 40×/0.65(S)	●
	Plan 100×/1.25(S)(Oil)	●
Trinocular Head	Inclined 45°, Light Distribution:20:80	●
Mechanical Stage	Travel Range: 50mm×40mm Size: 180mm×165mm	●
Stage Plate	1 (Φ10)	●
	2 (Φ20)	●
	3 (Φ40)	○
Set of Stage Plate		●
Main Body		●
Polarizing Slide		●
Analyzing Slide		●
Illuminator	6V/30W Halogen Lamp Spare Lamp: 2 pcs Spare Fuse: 2 pcs	●
Filters	Blue, Green, Grey, White	●
1× C-mountr	Available for Digital Camera	○
0.6× C-mount	Available for Vidicon	○
Available for Digital Camera		○
Micrometer	Division 0.01mm	○

Note: “●” in the table is standard attachment. “○” is optional accessories.

Design change: To keep pace with technological advances, we have reserved the right to make design modification and changes without notice.



Joy with your eyes

Chongqing MIC Technology Co., Ltd

Website: [www.micscope.com](http://www.micscope.com)

Email: [info@micscope.com](mailto:info@micscope.com)

TEL: +86-13436078184; FAX: +86-23-63913139