

## Y-ZM1 Microscope

	Excellent infinite Optical System. High Magnification Objective with Long Working Distance (up to 3.7mm for 40X Objective)
	Comfortable and Enjoyable Operation With low and Repetitive Position Control knob
	Triangle Structure Make it much steadier and less working Space Occupation
	Easy, Safe and Reliable Lamp Replacement Device

## Y-ZM1 Microscope

		Y-ZM1	
Viewing Head	Binocular Head	●	
	Trinocular Head	○	
Eyepiece	WF10×/18	●	
	WF8×/18	○	
	WF10×/18 with 0.1mm Micrometer & Cross Hair	○	
	WF12.5×/15	○	
Objective	Infinitive plan achromatic objectives	4×/0.1/∞/- WD25.4mm	●
		10×/0.25/∞/- WD11.0mm	●
		20×/0.40/∞/- WD6.0mm	○
		40×/0.60/∞/0 WD3.7mm	●
		40×/0.65/∞/- WD0.78mm	○
		50×/0.75/∞/- WD0.32mm	○
		80×/0.90/∞/- WD0.2mm	○
		100×/0.90/∞/0 WD0.7mm	○
Nosepiece	Quadruple Nosepiece	●	
Stage	Double layer mechanical stage 150×140mm/75×50mm	●	
Focusing	Coaxial Fine & Coarse Adjustment	●	
Illumination	Kohler Illumination, Epi-illuminator with iris aperture diaphragm and iris field diaphragm, 6V/20W Halogen Lamp, Brightness Adjustable	●	
	Diascopic Illumination, 6V/20W Halogen Lamp, Brightness Adjustable	○	
Filters	Blue, Green, Yellow & Ground Glass	●	
Polarization Set		○	
Specimen Presser	For metallurgical specimen preparation	○	

Note: “●” in the table is standard attachment. “○” is optional accessories.

Design change: To keep pace with technological advances, we have reserved the right to make design modification and changes without notice.



Joy with your eyes

Chongqing MIC Technology Co., Ltd

Website: [www.micscope.com](http://www.micscope.com)

Email: [info@micscope.com](mailto:info@micscope.com)

TEL: +86-13436078184; FAX: +86-23-63913139