

BS203ih Phase Contrast Microscope

Items Spec. and parameter	Model	BS203ih	Order No.
Eyepiece	WF10x/Ø18	●●	EP10g
	WF16x/Ø10	○○	EP16b
	4x Centering eyepiece	●	ECi
Plan Objectives	Infinite Plan Objectives 40x/0.1	●	B-OPi4a
	Infinite Plan Objectives 10x/0.25	●	B-OPi10a
	Infinite Plan Objectives 40x/0.65	●	B-OPi40a
	Infinite Plan Objectives 100x/1.25 (oil)	●	B-OPoi100a
	20x plan modulation contrast objective	●	B-OPih20
	40x plan modulation contrast objective	●	B-Opih40
Main body	Quadruple Nosepiece Coaxial Coarse and Fine Focusing Knobs (fine knob 0.002mm) Mechanical Stage (Size: 135×130mm, Travel Range: 50x75mm)	●	A14b
Binocular Head	Seidentopf Binocular Head Inclined 30°, Rotatable 360° Interpupillary Distance: 52-75mm	●	TBR20
Hofmann Modulation Contrast Condenser	NA=1.25	●	C32
Collector		●	COL03
Base	Bulb: Halogen Tungsten Lamp 6V20W, Fuse: Ø5X20/2A, Input voltage: AC85V-265V, 50HZ/60HZ, Output voltage: DC1.4V-6V Maximum output power 20W	●	B32

Note: “●” in the table is standard attachment. “○” is optional accessories.

Design change: To keep pace with technological advances, we have reserved the right to make design modification and changes without notice.



Chongqing MIC Technology Co., Ltd
Website: www.micscope.com
Email: info@micscope.com

TEL: +86-13436078184; FAX: +86-23-63913139