

IBE-1000 Microscope

Items	IBE1000	Order No.
10× WF eyepiece FN=20	●●	EW10
16× WF eyepiece FN=14	●●	EW16
10× plan reticle eyepiece FN=18	●	EPr10b
Centering eyepiece	●	ECi
2.5× long working distance objective NA=0.07; WD=28.89	○	B-OLPr2.5
4× long working distance objective NA=0.10 ;WD=16.50	○	B-OLPr4
10× long working distance objective NA=0.25; WD=7.90	●	B-OLPr10
25× long working distance objective NA=0.40; WD=5.00	●	B-OLPr25
40× long working distance objective NA=0.65; WD=3.00	●	B-OLPr40
10× long working distance positive phase contrast objective NA=0.25 ; WD=7.90	●	B-OLPpr10
25× long working distance positive phase contrast objective NA = 0.40 ;WD=5.00	●	B-OLPpr25
40× long working distance positive phase contrast objective NA = 0.65 ;WD=3.00	●	B-OLPpr40
10× long working distance negative phase contrast objective NA = 0.25 ;WD=7.90	○	B-OLPnr10
25× long working distance negative phase contrast objective NA=0.40 ; WD=5.00	○	B-OLPnr25
40× long working distance negative phase contrast objective NA = 0.65 ;WD=3.00	○	B-OLPnr40
long working distance phase contrast condenser NA=0.30;WD=55.00	●	C15
Mechanical stage: size:238×160; moving range: 50×50	●	S28
Sliding trinocular head, Inclined 45°, Interpupillary Distance: 55-75mm	●	18X.30
Photography tube group	○	MP04
Photography adaptor	○	TV04
Digital Photography device	○	CU43-2+DV02+camera
0.01mm micrometer	○	

Note: “●” in the table is standard attachment. “○” is optional accessories.

Design change: To keep pace with technological advances, we have reserved the right to make design modification and changes without notice.



Joy with your eyes

Chongqing MIC Technology Co., Ltd

Website: www.micscope.com

Email: info@micscope.com

TEL: +86-13436078184; FAX: +86-23-63913139