

## BA2000ip Microscope

	Order NO.	Item	BA2203ip
1	EW10	WF 10×Eyepiece FN=20	●●
2	EP16a	16×Plan Eyepiece FN=10	○○
3	EW16	WF 16× Eyepiece FN=14	○○
4	EPr10b	10×Reticule Eyepiece FN=18	●
5	B-OPi10a	10× Infinite Plan Achromatic Objective NA=0.25; WD=4.81	●
6	B-OPi20a	20× Infinite Plan Achromatic Objective, NA=0.40; WD=2.58	●
7	B-OPi40a	40× Infinite Plan Achromatic Objective, NA=0.66; WD=0.82	●
8	B-OPoi100a	100× Infinite Plan Achromatic Objective (oil), NA=1.25; WD=0.36	○
9	TTR11	Seidentopf Binocular Head, 30°inclined; interpupillary distance; 55~75	●
10	PA07	Middle Tube Group, Middle Tube Group: movable optical path, build in analyzer, 360° rotatable with Bertrand lens switch table, $\Lambda$ wave-plates, $\Lambda$ wave-plates and quartz fid ( I -IV )	●
11	SP01	Rotatable Stage, 360° rotatable with fixed thumb, lattice-value 1°, precision 6', size: $\phi$ 152	●
12	C28	Condenser NA=1.25	●
13		Nosepiece Quadruple, build in, center adjustable, inside inclined.	●
14	PO09	Polarizer 360° rotatable	●
15	PA06a.1	$\Lambda$ wave-plates	●
16	PA06a.2	$\Lambda$ wave-plates	●
17	PA06a.3	Quartz fid	●
18	A08.3	Base Tungsten lamp ; 6V/20W; wide voltage input	●
19		0.01 micrometer	○
20		Moving scale	○

- ★Infinite optical system with perfect optical characteristic
- ★Integrated frame design with stable operation quality
- ★Single polarizing observation, orthogonal polarization and conoscopic observation is available.
- ★Analysis and test for all kinds of mineral sample and other samples under polarizing.

Note: “●” in the table is standard attachment. “○” is optional accessories.

Design change: To keep pace with technological advances, we have reserved the right to make design modification and changes without notice.



Chongqing MIC Technology Co., Ltd  
 Website: [www.micscope.com](http://www.micscope.com)  
 Email: [info@micscope.com](mailto:info@micscope.com)  
 TEL: +86-13436078184; FAX: +86-23-63913139