

MIC-P3 Microscope

Specification		Model
		MIC-P3
Viewing head	Compensation Free Trinocular Head. Inclined at 30° (50mm -75mm)	●
Eyepiece	WF10 × /18mm	
	10 × gridding eyepiece	●
	10 × with reticule	●
	WF10 ×/ 22mm	●
Objective	Infinite plan-objective: 4 × ,10 × ,25 × , 40 × (S),63 × (S)	●
Stage	Round revolving stage	
	Stage dia: Φ 160mm	●
	Rotatable: 360°	
Condenser	Swing type N.A1.25 Abbe condenser with iris diaphragm & filter	●
Focusing	Coaxial coarse&fine focusing adjustment with rack and pinion mechanism Fine focusing scale Value 0.002mm	●
Collector	High brightness light collector	●
Light source	Halogen tungsten bulb 6V/30W.AC 85V-230V	●
	Adjustable brightness	●
Optional Accessory	Eyepieces: WF16x,WF20x,P16x 1.3Mega、 2.0Mega、 3.0Mega、 5.0Mega Pixels CMOS Digital camera eyepiece	

Note: “●” in the table is standard attachment. “○” is optional accessories.

Design change: To keep pace with technological advances, we have reserved the right to make design modification and changes without notice.



Chongqing MIC Technology Co., Ltd
Website: www.micscope.com
Email: info@micscope.com

TEL: +86-13436078184; FAX: +86-23-63913139