

MIC-S1 Stereo Microscope



MIC-S1-1A



MIC-S1-2B



MIC-S1-3C



MIC-S1-5C



MIC-S1-6C-W



MIC-S1-10



MIC-S1-20



MIC-S1-L

MIC-S1- series stereo microscopes are of fully coated excellent optical system providing high resolution and crisp erect image.

They are widely used in electric factories, schools laboratories, sculpture and families and so on.

		MIC-S1-1A	MIC-S1-1B	MIC-S1-1C
Head	Straight Binocular Viewing Head	●	●	●
Eyepiece	WF5×	○	○	○
	WF10×	●	●	●
	WF15×	○	○	○
	WF20×	○	○	○
Objective	Inserted Objective	1×	○	○
		2×	●	●
		3×	○	○
		4×	○	○
		6×	○	○
Working Distance	80mm	●	●	●
Illumination	No Light	●		
	Epi-illumination 12V/ 10W		●	●
	Diascopic Illumination 12V/ 10W			●



Joy with your eyes

Chongqing MIC Technology Co., Ltd

Website: www.micscope.com

Email: info@micscope.com

TEL: +86-13436078184; FAX: +86-23-63913139

MIC-S1 Stereo Microscope

		MIC-S1-2A	MIC-S1-2B	MIC-S1-2C	MIC-S1-3A	MIC-S1-3B	MIC-S1-3C	
Head	Inclined Binocular Viewing Head, Unrotatable	●	●	●	●	●	●	
Eyepiece	WF5×	○	○	○	○	○	○	
	WF10×	●	●	●	●	●	●	
	WF15×	○	○	○	○	○	○	
	WF20×	○	○	○	○	○	○	
Objective	Inserte d Objecti ve	1×	○	○	○			
		2×	●	●	●			
		3×	○	○	○			
		4×	○	○	○			
		6×	○	○	○			
	Objecti ve	2×、 4×				●	●	●
		1×、 2×				○	○	○
		1×、 3×				○	○	○
Working Distance	80mm	●	●	●				
	57mm				●	●	●	
Illuminat ion	No Light	●			●			
	Transmitted light 6V/ 12W		●	●		●	●	
	Incident light 6V/ 12W			●			●	

MIC-S1 Stereo Microscope

		MIC-S1-5A	MIC-S1-5B	MIC-S1-5C	MIC-S1-6C-W
Head	Inclined Binocular Head, 360° Rotatable	●	●	●	●
Eyepiece	WF5×	○	○	○	○
	WF10×	●	●	●	●
	WF15×	○	○	○	○
	WF20×	○	○	○	○
Objective	Inserted Objective	2×、4×	●	●	●
		1×、2×	○	○	○
	1×、3×	○	○	○	○
Working Distance	76mm	●	●	●	●
Illumination	No Light	●			
	Epi-illumination 12V/ 10W		●	●	●
	Diascopic Illumination 12V/ 10W			●	●

		MIC-S1-10	MIC-S1-20	
Head	Inclined Binocular Head	●	●	
Eyepiece	WF5×	○	○	
	WF10×	●	●	
	WF15×	○	○	
	WF20×	○	○	
Objective	Objective	2× & 4×	●	●
		1× & 2×	○	○
		1× & 3×	○	○
	Auxiliary Objective	0.5×	○	○
		0.75×	○	○
		1.5×	○	○
	2×	○	○	
Working Distance	72mm	●	●	
Illumination	No Light	●		
	Epi-illumination 6V/ 12W		●	
	Diascopic Illumination 6V/ 12W		●	

MIC-S1 Stereo Microscope

		MIC-S1-L
Head	Inclined Binocular Head	
Eyepiece	WF5×	
	WF10×	
	WF15×	
	WF20×	
Objective	Inserted Objective	1×
		0.5×
Working Distance	230mm (Inserted Objective 1×)	
	250mm (Inserted Objective 0.5×)	
Illumination	Incandescent Lamp 12V/ 10W	

Note: “●” in the table is standard attachment. “○” is optional accessories.

Design change: To keep pace with technological advances, we have reserved the right to make design modification and changes without notice.



Chongqing MIC Technology Co., Ltd
 Website: www.micscope.com
 Email: info@micscope.com

TEL: +86-13436078184; FAX: +86-23-63913139