

MIC-S13 Stereo Microscope

MIC-S13 is the new large viewing field microscope. It produces sharp and stereoscopic images. It is to be operated easily and comfortably. Widely applied to biological dissection, electronics, minerals.

Features

- Imaging from optical entity; providing with various objectives of 2×、4×、6×、8×、10×、15× magnification. The objectives can be conveniently changed; owning humanistic operating system.
- For the long time viewing without visual weariness, the freely and comfortably operating can be beyond the visual limit, and improve the working efficiency.
- Possessing the super-long working distance super-large field of viewing and depth of field, so that the operator hasn't any limited.

Application to electronic industry precision project plastic shop medical diagnosis and biological medicine public security system and so on.

Specification & Objective

Magnification	Working distance(mm)	Viewing field for object space(mm)	Equipment
2×	208	68	●
4×	98	34	●
6×	80	22.7	●
8×	58	17	○
10×	46	13.6	○
15×	50	9.1	○
6× (SL)	115	22.7	○



Floor model wide-field
Microscope

"●" is standard equipment. "○" is optional equipment.

Height: 1800mm

Free Viewing Body, Strength Range: 850mm,
easily observe large area objective.

Stable stand, free sliding, long vertical distance is
580mm with boost up and down.



Chongqing MIC Technology Co., Ltd
Website: www.micscope.com
Email: info@micscope.com

TEL: +86-13436078184; FAX: +86-23-63913139

MIC-S13 Stereo Microscope

Universal block with tension adjustment in bottom,
it can free to move or fix so that to observe in a
good position.

Application range

- Electronics: PCB Setting, QC, Screening, Rework & Weld
- Precision project & Plastic: QC, Micro-welding, Micro-machining Injecting & Moulding
- Medical & Dentistry Equipment
Manufacturing: Finish machining Setting, Fine Dressing, Color Matching, Rework & Overhaul
- Biological medicine: Sample-Making, dissecting, Dyeing, Microscopy
- Public Security system



Optional LED Illumination

MIC-S13 Stereo Microscope



Note: “●” in the table is standard attachment. “○” is optional accessories.

Design change: To keep pace with technological advances, we have reserved the right to make design modification and changes without notice.

CE ISO9001
ISO14001

MIC[®]
Joy with your eyes

Chongqing MIC Technology Co., Ltd

Website: www.micscope.com

Email: info@micscope.com

TEL: +86-13436078184; FAX: +86-23-63913139