

## MIC-S3 Stereo Microscope

		MIC-S3
Drawing Tube	Binocular Viewing Head, Inclined at 45°	●
	Trinocular Viewing Head, Inclined at 45°	○
Eyepiece	Extra wide field eyepiece EW10×/22	●
	WF15×	○
	WF20×	○
Zoom Objective	0.8×-5×	●
Zoom Ratio	6.3: 1	●
Working Distance	115mm	●
Interpupillary Distance	57-75mm	●
Auxiliary objective	0.5×	○
	2×	○
Glass Insert	Glass Insert, Diameter 100mm	●
Illumination	Incident Illumination 100V-240V/ LED	●
	Transmitted Illumination 100V-240V/ LED	●
Photo Attachment		○
Video Attachment		○

Note: “●” in the table is standard attachment. “○” is optional accessories.

Design change: To keep pace with technological advances, we have reserved the right to make design modification and changes without notice.



Chongqing MIC Technology Co., Ltd  
 Website: [www.micscope.com](http://www.micscope.com)  
 Email: [info@micscope.com](mailto:info@micscope.com)

TEL: +86-13436078184; FAX: +86-23-63913139