

MIC-S7 Stereo Microscope

| Specification | | Model | | |
|--------------------|--|---------|---------|---------|
| | | MIC-S7A | MIC-S7B | MIC-S7C |
| Viewing Head | Binocular Inclined at 45°, 360° Rotatable (55mm-75mm) | ● | ● | ● |
| Eyepiece | WF10X/20mm | ● | ● | ● |
| Objective | 1X、2X Working distance:100mm | ● | | |
| | 1X、3X Working distance:100mm | | ● | |
| | 2X、4X Working distance:100mm | | | ● |
| Illumination | Incident Light Halogen lamp 12V10W or 12V15W; Transmitted Light Halogen lamp; 12V10W or 12V15W | ● | ● | ● |
| Optional Accessory | Eyepiece: SWF15X/15 SWF20X/11 Ring Fluorescent Light; 1.3Mega pixels Digital camera eyepiece; LED illumination | | | |

Note: “●” in the table is standard attachment. “○” is optional accessories.

Design change: To keep pace with technological advances, we have reserved the right to make design modification and changes without notice.



Chongqing MIC Technology Co., Ltd
 Website: www.micscope.com
 Email: info@micscope.com

TEL: +86-13436078184; FAX: +86-23-63913139