

SD200 Stereo Microscope

Order No.	Item	SD101	SD102	SD111	SD112
ES10	Wide Field Plan Eyepiece WF10X/22mm	●	●	●	●
ES15	Plan Eyepiece WF15X/16mm	●	●	●	●
ES20	Plan Eyepiece WF20X/13mm	●	●	●	●
ES25	Plan Eyepiece WF25X9mm	○○	○○	●●	●●
ES30	Plan Eyepiece WF30X/7mm	●	●	●	●
8T.41X	Plan Objective W.D.134mm	○○	○○	●●	●●
8T.52X	Plan Objective W.D.134mm	●●	●●	●●	●●
8T.63X	Plan Objective W.D.110mm	○○	○○	○○	○○
8T.74X	Plan Objective W.D.110mm	○○	○○	○○	○○
SD200	Head (inclined 30°)	○○	○○	○○	○○
	Main Body (Mount SizeΦ80mm)	○○	○○	○○	○○
	Halogen Tungsten Reflecting Lamp: 6V15W (Inclined Illumination) Input: AC220V/50HZ	●		●	
	Fluorescent Lamp 220V7W (Transmission Illumination)	●	●		
	Halogen Tungsten Lamp 6V15W (Transmission Illumination)	○	○		
TV11	CCD Photo Eyepiece	●	●	●	●
	CCD Video Camera	○	○	○	○
	Ring Fluorescent Lamp	○	○	○	○
	LED Ring Fluorescent Lamp	○	○	○	○
	Black And White Plastic Stage Size: 250×220mm/Ø95	○	●	○	●
	Ground Glass Stage Size: 250×220mm/Ø95	●	○	●	○

Note: “●” in the table is standard attachment. “○” is optional accessories.

Design change: To keep pace with technological advances, we have reserved the right to make design modification and changes without notice.



Joy with your eyes

Chongqing MIC Technology Co., Ltd

Website: www.micscope.com

Email: info@micscope.com

TEL: +86-13436078184; FAX: +86-23-63913139