

## XTD-501 Stereo Microscope

Name	Stereo Microscope
Model	XTD-501
Specification	1. Magnification: 20×--40× 2. Eyepiece: WF10×/ F.N 20 3. Objective: 2×/ W.D 93, 4×/ W.D 66 4. Rotatable Binocular (360°) 5. Illumination: Transmission and Reflect Light
Optional	Objective 1×/ W.D 80 Objective 3×/ W.D 69 Eyepiece WF 15×/ F.N 13 Eyepiece WF 20×/ F.N 10
Characteristic and Function	1. This microscope can produce clear image with balanced illumination light. It is easy to operate with comfortable observation, delicate construction and smart shape. 2. This microscope can be used in the inspection, assembly and repairing of the electron components and slim precise accessory in electron industry, apparatus and instrument, precise mechanism industry. It can also be used in farming, spinning, medical dissection, geology, archaeology and carving for inspection, observation and even cultural relic identification.

Note: “●” in the table is standard attachment. “○” is optional accessories.

Design change: To keep pace with technological advances, we have reserved the right to make design modification and changes without notice.



Chongqing MIC Technology Co., Ltd  
 Website: [www.micscope.com](http://www.micscope.com)  
 Email: [info@micscope.com](mailto:info@micscope.com)

TEL: +86-13436078184; FAX: +86-23-63913139