

XTL-3 Stereo Microscope

Item		Model	Configuration			Order NO.
			XTL-3	XTL-3A	XTL-3B	
Main Body	XTL-3		●			
	XTL-3A			●		
	XTL-3B				●	
Objective	Zoom objective: 0.7X-3X(zoom ratio: 4.28:1)		●	●	●	
Eyepiece	Kellner 10X-18		●●	●●	●●	EK10X
	Kellner15X-15		○○	○○	○○	EK15X
	Kellner10X (reticle)		○	○	○	EK10 (K)
Lamp	Halogen Bulb 6V 10W					
Stage (ø 85)						3T.1-5A
Eye cap			○○	○○	○○	2T-2

Note: “●” in the table is standard attachment. “○” is optional accessories.

Design change: To keep pace with technological advances, we have reserved the right to make design modification and changes without notice.

CE **ISO9001**
ISO14001



Chongqing MIC Technology Co., Ltd
Website: www.micscope.com
Email: info@micscope.com

TEL: +86-13436078184; FAX: +86-23-63913139