

# XTL-3B Stereo Microscope

## XTL-3B Stereo Microscope Standard Outfits

Item	XTL-3B
Objective 1x~4x; WD=94	●
Wide Field Eyepiece WF10X/22mm	●●
Wide Field Eyepiece WF20X/14mm	○○
Video Eyepiece	○
Seidentopf Trinocular Head / Inclined 45° Interpupillary Distance: 55-75mm	●
Stage, Glass	●
Clip	●●
Project light Halogen Tungsten Lamp 6V10W	●
Transmit light Halogen Tungsten Lamp 6V10W	●

Note: “●” in the table is standard attachment. “○” is optional accessories.

Design change: To keep pace with technological advances, we have reserved the right to make design modification and changes without notice.



Chongqing MIC Technology Co., Ltd

Website: [www.micscope.com](http://www.micscope.com)

Email: [info@micscope.com](mailto:info@micscope.com)

TEL: +86-13436078184; FAX: +86-23-63913139