

ZSA300 Stereo Microscope

ZSA300 Series Stereo Microscope Standard Outfits
(Zoom Range: 0.7~5X/WD129.5mm)

Order No.	Item	ZSA301	ZSA302
8T. 1a	Main Body / Mount Dia. Ø80mm	●	●
18T	Main Body / /Mount Dia. Ø80mm		
ES10	Wide Field Eyepiece WF10X/22mm	●●	●●
ES15	Wide Field Eyepiece WF15X/16mm	○○	○○
ES20	Wide Field Eyepiece WF20X/13mm	○○	○○
ES25	Wide Field Eyepiece WF25X/9mm	○○	○○
ES30	Wide Field Eyepiece WF30X/7mm	○○	○○
18T.4	0.5x Big Objective	○	○
18T.5	1.5x Big Objective	○	○
18T.6	2x Big Objective	○	○
18T	Binocular Head / Inclined 30°	●	
18TA	Trinocular Head / Inclined 30°		●
18TB	Binocular Head / Inclined 30°		
	White and Black Plastic Stage / Ø95mm	○	○
	Ground Glass Stage / Ø95mm	●	●
TV12	CCD TV Attachment		●
	Halogen Tungsten Lamp 6V15W (Inclined Illumination)	●	●
	Fluorescent Lamp 220V7W (Transmission Illumination)	●	●
	Halogen Tungsten Lamp 6V20W (Transmission Illumination)	○	○
	Halogen Tungsten Lamp 6V30W(Ribbon with Coaxial illumination)		
	Input: AC220V/50Hz Fuse Tube Ø5×20 2A		

Note: “●” in the table is standard attachment. “○” is optional accessories.

Design change: To keep pace with technological advances, we have reserved the right to make design modification and changes without notice.



Chongqing MIC Technology Co., Ltd
Website: www.micscope.com
Email: info@micscope.com

TEL: +86-13436078184; FAX: +86-23-63913139